

Zipper® Leg Sleeve Expander Zipper® Arm Sleeve Expander

Zipper® Arm/Leg Sleeve expander is a fabric strip, that has two Zipper fasteners on each side along the strip




Use the Zipper® Arm/Leg Sleeve expander with Mego Afek Garments:
Zipper® Leg Sleeve | Zipper® arm Sleeve

Instructions for Wearing:

1. Open the zipper fastener and wrap over your leg or over arm
2. Close the zipper fasteners, closing the sleeve.
3. **If the sleeve is still too small to cover the limb properly, open the sleeve zippers completely and attach an Expander to the Sleeve's Zippers.**

Zipper® Leg Sleeve Expander Zipper® Arm Sleeve Expander

Availability Chart:



Cat. No. (Expander)	Fits Sleeve	Zipper Length A	Circumference increase
L4159370	ZIP 3-60, Pedi 4-55	50	6.5
L4159290	ZIP 3-70	62	6.5
L4150680	ZIP 1-75S, ZIP 1-75M, ZIP 1-75L, ZIP 1-80S, ZIP 1-80M, ZIP 1-80L, ZIP 2-75, ZIP 2-80, ZIP 3-75	68	6.5
L4150750	ZIP 2-85, ZIP 3-85, ZIP 4-85	75	6.5
L4150850	ZIP 3-90	85	6.5

1. All Dimensions are in cm (Centimeters)
2. General tolerance: ± 1 cm